

TEST INKS BLUE AND TEST PENS BLUE

from 18.4 to 105 mN/m / toxic from 24 to 56 mN/m

Blue coloured Test Inks are special testing liquids in ranges of defined surface tension (in accordance with ISO 8296 / ASTM 2578 / DIN 53364). The surface tension of a substrate is checked by simply applying the Test Ink to the surface. The blue Test Inks are available in a range from 18.4 to 105 mN/m. Test Pens filled with blue Test Ink are available from 28 - 70 mN/m.



TEST PENS

- from 28 60 mN/m in steps of one by one from 62 70 mN/m in steps of two by two
- with accuracy +/- 1.0 mN/m
- observation time is 2 Sec.
- simple handling
- no spilling

TEST INKS in bottles

- available from 18.4 to 105 mN/m (18.4, 76, 84, 90, 105mN/m colourless)
- with accuracy +/- 0.5 mN/m
- observation time is 2 Sec.
- bottles containing 10, 100, 250 und 1000 ml or in sets of 7 bottles of 10 ml each

Germany