

arcotest[®] ORGANIC Test Inks and arcotest[®] ORGANIC Test Pens

available from 31 to 46 mN/m

arcotest[®] ORGANIC Test Inks are special testing liquids in ranges of defined surface tension – in green colour. They were developed to obtain non-toxic test inks that do not require labelling. The application and handling of arcotest[®] ORGANIC is based on DIN 53364 re. ISO 8296. The ORGANIC Test Inks are available in bottles or in pens from 31 to 46 mN/m (Dyn/cm).



arcotest[®] ORGANIC in bottles

- from 31 to 46 mN/m
- accuracy +/- 0.5 mN/m
- observation time is 2 sec.
- bottles containing 10, 100 and 250 ml or in sets of 7 bottles of each 10 ml

arcotest[®] ORGANIC Test Pens

- \bullet from 31 to 46 mN/m
- accuracy +/- 1.0 mN/m
- observation time is 2 sec.
- simple handling
- no spillage
- line width 5 mm
- pen sets available in sets of 4, 6 and 8 pcs. or as single pens, 5 ml

Not poisonous Not harmful to the environment Not subject to labelling Not harmful to health

arcotest[®] ORGANIC Test Pens Jumbo

- 32 to 46 mN/m in steps of 2
- accuracy +/- 1.0 mN/m
- observation time is 2 sec.
- simple handling
- no spilling
- line width 15 mm
- optimal for testing large film areas
- pen sets available in sets of 4, 6 and 8 pcs.
- or as single pens, 17 ml

Managing Director: Fritz P. Bloss Amtsgericht Mannheim HRB 511684 VAT DE 813775094 **arcotest GmbH** Rotweg 25 71297 Moensheim Germany www.arcotest.info E-mail: info@arcotest.info Phone +49(0)7044 - 902 270 Fax +49(0)7044 - 902 269

