

Test Inks ORGANIC in bottles

from 30 to 46 mN/m // with accuracy measurement ± 0.5 mN/m // without brush // not subject to labelling



STANDARD set

32 34 36 38 40 42 44 mN/m
7 glass bottles of 10 ml
art. no. 40.20000.4 without brush



OPTIONAL set of 7

30 to 46 mN/m
7 glass bottles of 10 ml of your choice
art. no. 40.20001.4 without brush



10 ml

30 to 46 mN/m
standard and special
art. no. 40.201XX.4 without brush



PET 10 ml

30 to 46 mN/m
Plastic bottle, especially for use in the food sector
art. no. 40.204XX.4 without brush



100 ml

30 to 46 mN/m
standard and special
art. no. 40.202XX.0
art. no. 40.202XX.0 (SE)



250 ml

30 to 46 mN/m
standard and special
art. no. 40.203XX.0
art. no. 40.203XX.0 (SE)



arcoweb®

38 mN/m
Application width 25 mm
art. no. 40.20438.0



arcospray 38®

38 mN/m
labeling required
for small and large areas
art. no. 40.80238.0

Test Pens ORGANIC

from 30 to 46 mN/m // with accuracy measurement ± 1.0 mN/m // not subject to labelling



STANDARD set of 8

30 32 34 36 38 40 42 44 mN/m
8 test pens
art. no. 40.25000.0



OPTIONAL set of 8

30 to 46 mN/m
8 test pens of your choice
art. no. 40.25001.0



OPTIONAL set of 6

30 to 46 mN/m
6 test pens of your choice
art. no. 40.25002.0



OPTIONAL set of 4

30 to 46 mN/m
4 test pens of your choice
art. no. 40.25003.0



TEST PEN

30 to 46 mN/m
standard and special
art. no. 40.251XX.0

Test Pens ORGANIC Jumbo

from 30 to 46 mN/m in steps of 2 // with accuracy measurement ± 1.0 mN/m // not subject to labelling



STANDARD set of 8

30 32 34 36 38 40 42 44 mN/m
8 test pens
art. no. 40.26000.0



OPTIONAL set of 8

30 to 46 mN/m
8 test pens of your choice
art. no. 40.26001.0



OPTIONAL set of 6

30 to 46 mN/m
6 test pens of your choice
art. no. 40.26002.0



OPTIONAL set of 4

30 to 46 mN/m
4 test pens of your choice
art. no. 40.26003.0



TEST PEN Jumbo

30 to 46 mN/m
art. no. 40.261XX.0

Test Inks PINK in bottles

from 22 to 60 mN/m // non-toxic
with accuracy measurement ± 0.5 mN/m

STANDARD: 30 32 34 36 38 40 42 44 46 48 50 52 54 56 58 60 mN/m
SPECIAL(SE): 22 24 26 28 31 33 35 37 39 41 43 45 47 49 51 53 55 57 59 mN/m



STANDARD set A

32 34 36 38 40 42 44 mN/m

7 glass bottles of 10 ml
art. no. 40.60000.0 with brush
art. no. 40.60000.4 without brush



OPTIONAL set of 7

22 to 60 mN/m

7 glass bottles of 10 ml
art. no. 40.60001.0 with brush
art. no. 40.60001.4 without brush



10 ml

22 to 60 mN/m

standard and special
art. no. 40.601XX.0 with brush
art. no. 40.601XX.4 without brush



PET 10 ml

32 to 44 mN/m

Plastic bottle, especially for use in the food sector
art. no. 40.700XX.0 with brush
art. no. 40.700XX.4 without brush



100 ml

22 to 60 mN/m

standard and special
art. no. 40.602XX.0
art. no. 40.602XX.0 (SE)



250 ml

22 to 60 mN/m

standard and special
art. no. 40.603XX.0
art. no. 40.603XX.0 (SE)

Test Pens PINK

from 22 to 60 mN/m // non-toxic
with accuracy measurement ± 1.0 mN/m

STANDARD: 30 32 34 36 38 40 42 44 mN/m
SPECIAL (SE): 22 24 26 28 31 33 35 37 39 41 43 45 - 60 mN/m



STANDARD set of 8

30 32 34 36 38 40 42 44 mN/m

8 test pens
art. no. 40.45001.0



OPTIONAL set of 8

22 to 60 mN/m

8 test pens of your choice
art. no. 40.45000.0



OPTIONAL set of 6

22 to 60 mN/m

6 test pens of your choice
art. no. 40.45002.0



OPTIONAL set of 4

22 to 60 mN/m

4 test pens of your choice
art. no. 40.45003.0



TEST PEN

22 to 60 mN/m

standard and special
art. no. 40.451XX.0

Test Pens PINK Jumbo

from 30 to 50 mN/m // non-toxic // with accuracy measurement ± 1.0 mN/m



STANDARD set of 8

30 32 34 36 38 40 42 44 mN/m

8 test pens
art. no. 40.46000.0



OPTIONAL set of 8

30 to 50 mN/m

8 test pens of your choice
art. no. 40.46001.0



OPTIONAL set of 6

30 to 50 mN/m

6 test pens of your choice
art. no. 40.46002.0



OPTIONAL set of 4

30 to 50 mN/m

4 test pens of your choice
art. no. 40.46003.0



TEST PEN Jumbo

30 to 50 mN/m

art. no. 40.461XX.0

Test Inks BLUE in bottles

from 18 to 105 mN/m // toxic from 24 to 57 mN/m // with accuracy measurement ± 0.5 mN/m // 30 to 72 mN/m according to ISO 8296 (DIN 53364 and ASTM)

STANDARD: 28 32 35 38 41 44 48 56 mN/m
 SPECIAL (SE): 18 (colourless) 20 22 24 26 29 30 31 33 34 36 37 39 40
 42 43 45 46 47 49 50 51 52 53 54 55 57 58 60 62 64
 66 68 70 72 mN/m
 colourless: 74 76 84 90 105 mN/m



Set A

28 35 38 41 44 48 56 mN/m
 7 glass bottles of 10 ml
 art. no. 40.30001.0 with brush
 art. no. 40.30001.4 without brush



Set B

28 32 35 38 41 44 48 mN/m
 7 glass bottles of 10 ml
 art. no. 40.30000.0 with brush
 art. no. 40.30000.4 without brush



Set C

30 32 34 36 38 40 42 mN/m
 7 glass bottles of 10 ml
 art. no. 40.30003.0 with brush
 art. no. 40.30003.4 without brush



OPTIONAL set

18 (colourless) 20 to 72 mN/m
 7 glass bottles of 10 ml of your choice
 art. no. 40.30002.0 with brush
 art. no. 40.30002.4 without brush



10 ml

18 (colourless) to 72 mN/m
 standard and special
 art. no. 40.301XX.0 with brush
 art. no. 40.301XX.4 without brush



100 ml

18 (colourless) to 72 mN/m
 standard and special
 art. no. 40.302XX.0
 art. no. 40.302XX.0 (SE)



250 ml

18 (colourless) to 72 mN/m
 standard and special
 art. no. 40.303XX.0
 art. no. 40.303XX.0 (SE)



COLOURLESS 10 ml

74 76 84 90 105 mN/m
 special
 art. no. 40.301XX.0 (SE) with brush

Test PENS BLUE

from 28 - 60 mN/m in steps of one by one // from 62 - 72 mN/m in steps of two by two
 toxic from 24 to 57 mN/m // with accuracy measurement ± 1.0 mN/m



OPTIONAL set of 8

28 to 72 mN/m
 8 test pens of your choice
 art. no. 40.35001.0



OPTIONAL set of 6

28 to 72 mN/m
 6 test pens of your choice
 art. no. 40.35002.0



OPTIONAL set of 4

28 to 72 mN/m
 4 test pens of your choice
 art. no. 40.35003.0



TEST PEN

28 to 72 mN/m
 art. no. 40.351XX.0

QUICKTEST 38®



QUICKTEST 38®

38 mN/m approx.
quick check for polyolefins,
5 ml
art. no. 40.55100.0

RAPIDTEST 38®



QUICKTEST 38® Jumbo

38 mN/m approx.
quick check for polyolefins,
15 ml, 15 mm line width
art. no. 40.55100.4



RAPIDTEST 38®

38 mN/m approx.
quick check for polyolefins,
5 ml
art. no. 40.66100.0



RAPIDTEST 38® Jumbo

38 mN/m approx.
quick check for polyolefins,
15 ml, 15 mm line width
art. no. 40.66100.4

Case SET



arcotest® PINK case SET

Case content:
• PINK set A 32 – 44 mN/m
without brush
• arcotestCLEANER 250 ml
• cotton tipped applicators
art. no. 40.60000.8
art. no. 40.60001.8 (optional set)



arcotest® BLUE case SET

Case content:
• BLUE set A 28 – 56 mN/m
without brush
• arcotestCLEANER 250 ml
• cotton tipped applicators
art. no. 40.30001.8
art. no. 40.30002.8 (optional set)



arcotest® ORGANIC case SET

Case content:
• ORGANIC set 32-44 mN/m
without brush
• cotton tipped applicators
art. no. 40.20000.8
art. no. 40.20001.8 (optional set)

Case



TRANSPORT CASE, large

for 24 test ink bottles of 10 ml
without contents
art. no. 40.31900.0



TRANSPORT CASE, small

for 7 test ink bottles of 10 ml, with
storage for test pen (set of 4) or cotton
tipped applicators
without contents
art. no.40.31800.0 (blue)
art. no. 40.31800.5 (green)



arcotest® Plasma-Set

Case content:
• handheld plasma device
piezobrush® PZ3 for manual
surface treatment
• arcotest® 4 test pens of your choice
art. no. 40.00000.8

Accessories



COTTON-TIPPED APPLICATORS

approved, 100 piece/pack,
15 cm long
art. no. 40.31700.0



TEST LIGHT

for use in case of low contrast between
measuring surface and test ink
art. no. 40.31600.0



arcotestCLEANER

optimizes cleanliness
content 250 ml
art. no. 40.32000.0