

Test Inks PINK and Test Pens PINK

available from 22 to 60 mN/m

PINK coloured Test Inks are special testing liquids in ranges of defined surface tension. They were developed to get "non toxic" Test Inks . These PINK Test Inks are meant to substitute the BLUE coloured inks as they were declared to be toxic according to DIN 53364 / ISO 8296. The PINK Test Inks are available in bottles or in pens.



Test Inks PINK in bottles

- from 22 to 60 mN/m
- accuracy +/- 0,5 mN/m
- observation time 22 to 26 mN/m is 2 sec.
- observation time 28 to 44 mN/m is 4 sec.
- observation time 45 to 60 mN/m is 2 sec.
- bottles containing 10 or 100 or 250 mL
 or in sets of 7 bottles of each 10 mL

Test Pens PINK

- from 22 to 60 mN/m
- accuracy +/- 1,0 mN/m $\,$
- observation time 22 to 26 mN/m is 2 sec.
- observation time 28 to 44 mN/m is 4 sec.
- observation time 45 to 60 mN/m is 2 sec.
- simple handling
- · no spilling
- pen sets available in sets of 4, 6 and 8 pcs.
 or as single pens, 5 mL

Test Pens PINK Jumbo

- from 30 to 50 mN/m $\,$
- accuracy +/- 1,0 mN/m
- observation time 30 to 44 mN/m is 4 sec.
- observation time 46 to 50 mN/m is 2 sec.
- simple handling
- no spilling
- line width 15 mm
- optimal for testing large film areas
- pen sets available in sets of 4, 6 and 8 pcs.
 or as single pens, 17 mL

Managing Director: Fritz P. Bloss Amtsgericht Mannheim HRB 511684 IdNr. DE 813775094



Rotweg 25 71297 Moensheim Germany www.arcotest.info E-mail: info@arcotest.info Phone +49(0)7044 - 902 270 Fax +49(0)7044 - 902 269

