

arcotest®

Made in Germany

arcotest® ORGANIC Test Inks

available from 30 to 60 mN/m

arcotest® ORGANIC Test Inks are special testing liquids in ranges of defined surface tension – in green colour. They were developed to obtain non-toxic test inks that do not require labelling. The application and handling of arcotest® ORGANIC is based on DIN 53364 re. ISO 8296. The ORGANIC Test Inks are available in bottles, in pens or arcoweb® ORGANIC disposable cloths.



Not poisonous
Not harmful to the environment
Not subject to labelling
Not harmful to health

arcotest® ORGANIC in bottles

30 to 60 mN/m

- accuracy: ± 0.5 mN/m
- observation time: 2 sec.
- Available:
 - 10 ml bottles with or without brush
 - sets of 7×10 ml
 - 100 and 250 ml bottles, ink sampling with arcotest® cotton swabs

Test Pens arcotest® ORGANIC

30 to 60 mN/m

- accuracy: ± 1.0 mN/m
- observation time: 2 sec.
- line width: 5 mm
- simple handling
- no spilling
- Available:
 - sets of 4, 6 and 8 pcs.
 - single pen, 5 ml

arcoweb® Disposable Cloth

38 mN/m

- accuracy: ± 1.0 mN/m
- observation time: 2 sec.
- application width: 25 mm
- consistently small application volume
- optimal for testing large areas
- Available:
 - individual packages

Test Pens arcotest® ORGANIC Jumbo

30 to 46 mN/m in steps of 2

- accuracy: ± 1.0 mN/m
- observation time: 2 sec.
- line width: 15 mm
- simple handling
- optimal for testing large areas
- Available:
 - sets of 4, 6 and 8 pcs.
 - single pen, 5 ml

Managing directors:
Sebastian Bloss, Fritz P. Bloss
Amtsgericht Mannheim HRB 511684
VAT.No. DE 813775094

arcotest GmbH
Rotweg 25
71297 Moensheim
Germany

www.arcotest.info
E-mail: info@arcotest.info
Phone: +49(0)7044 - 902 270
Fax: +49(0)7044 - 902 269

